Prisma
APHTM

Prisma
APHTM

Manter fora do alcance de crianças
Mantener fuera del alcance de los niños
Keep out of reach of children

Uso profesional
Uso profesional
Professional use

Fabricado por/Ficheiro por/Manufactured by:
DENTSPRY Ind. e Com. Ltda. - Rua Alice Heron, 86 -
Petrópolis/RJ
CEP: 25665-010 - CNPJ 31.116.239/0001-55 -
www.dentsply.com.br

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USO PROFESIONAL

PRISMA LIGHT CURED APH HYBRID RESIN

Light Cured Composite Resin

PRESENTATION

REFILL with 01 Prisma APH Syringe 4g, in one of the colors: XL, L, LG, U-(LYG), LY, DG, DY, GB, A2, A3, A3.5, B3, D3, O-A2 O-A3.5

COMPOSITION

Aluminum borosilicate, pigments, Aerosil R-972, urethane-modified BIS-GMA resin, BHT, Camphorquinone, FluBlau Concentrate, Uvinul M40 and N-methyl diethanolamine.

INDICATIONS FOR USE

Prisma APH is a radiopaque restorative material for restoration of anterior and posterior teeth, activated by visible light. It should be used in conjunction with the Prime & Bond 2.1, XP Bond Adhesive Systems or other universal adhesive system. Prisma APH comes in syringes with 4 grams and 15 color options to choose from.

PRECAUTIONS/WARNINGS

- The Tooth Conditioner Gel contains phosphoric acid, which can cause burns. Avoid contact with oral tissues, skin and eyes. If accidental contact occurs, wash the area with a generous amount of water. For eyes, flush with plenty of water and seek medical assistance.
- The Dental Gel Conditioner should be dispensed easily: NEVER USE EXCESSIVE FORCE. Before using the Dental Gel Conditioner, make sure that the applicator needles are properly fixed to the syringe body. Test the flow of material away from the operative field.
- The applicator needles may be slightly curved to facilitate access to the cavities, however, avoid sharp angles since these can obliterate the needle cannula. After each use, remove the needle and replace the protective cover.
- Prisma APH and Prime & Bond 2.1 contain polymerizable monomers, which may cause skin sensitization (allergic contact dermatitis) in susceptible individuals. After contact, wash the area with soap and water. If sensitization or any other allergic reaction occurs, discontinue use.
- Prime & Bond 2.1 contains acetone, which is highly flammable. Keep away from open flames, do not smoke, do not inhale vapors.
- Volatile - Replace covers of the Prime & Bond 2.1 immediately after each use.
- Keep Prisma APH and Prime & Bond 2.1 away from sunlight and in a ventilated spot.
- Do not store at temperatures higher than 30° C.

STEP-BY-STEP INSTRUCTIONS

Correlation of colors between the Prism and Standard Scales:

Prisma Scale	Standard Scale	Prisma Scale	Standard Scale
XL	A1	LG	C1
- - -	A2	U (LYG)	C2
- - -	A3	DG	C3
- - -	A3,5	GB	C4
DY	A4	- - -	D3
L	B1	- - -	O-A2
LY	B2	- - -	O-A3,5
- - -	B3		

COLOR SELECTION

The color selection should be done after prophylaxis with Odahcam Prophylaxis Paste or Pumice and while the tooth is wet.

PREPARATION OF THE CAVITY

- 1 - Restorations of anterior teeth
Use conventional cavity preparations for Class III, IV and V cavities
- 2 - Restorations of posterior teeth
Use conventional cavity preparations, however, make a bevel at a cavosurface angle to increase the area of enamel etching. Prevent amalgam residue or other base material from going inside the cavity, as these can damage the correct polymerization of light curing restorative materials.

PULP PROTECTION (before etching)

For deep cavities whose floor has great proximity to the dental pulp, the application of a Dycal or Hydro C type calcium hydroxide liner is recommended.

PLACEMENT OF MATRIX

To anatomically restore anterior teeth and minimize any excess material common or preformed polyester matrices should be used. In posterior teeth, one should use metal or polyester matrices made with a cast, which allows you to obtain adequate proximal contacts.

ACID ETCHING

The prepared enamel must be clean and dry. Apply the Tooth Conditioner Gel to the enamel walls for 30 seconds and to the and dentin for a maximum of 15 seconds, rinse the cavity and dry with quick bursts of air, avoiding desiccation of the exposed dentin. The properly conditioned enamel will have a whitish matte look.

Specifications - chemical composition of the product:

Components			
Urethane-Modified Bis-GMA Resin	Pyrolytic Silanized Silica	N. Methyl Diethanolamine	Mineral Dyes
Aluminum borosilicate and Silanized Barium	Camphorquinone	Butylated hydroxytoluene	Colloidal Silica

APPLICATION OF PRIME & BOND 2.1

- 1 - Dispense 2-3 drops of Prime & Bond 2.1 in a clean dappen pot. Hermetically close the flask to prevent evaporation of the product.
- 2 - Apply a generous amount of Prime & Bond 2.1, allowing all cavity surfaces to be properly wet. The cavity surfaces must remain at rest for 20 seconds and may require additional applications of Prime & Bond 2.1
- 3 - Remove the excess bonding agent with a brief blast of air (5 seconds max). The surface should have a uniform brightness. If it does not, repeat the steps in items 1 and 3.
- 4 - Light cure Prime & Bond 2.1 for 10 seconds.
- 5 - Apply a second layer of Prime & Bond 2.1 over the entire cavity and immediately remove the excess with a blast of air.
Light cure for 10 seconds and proceed with the restoration.

APPLICATION OF PRISMA APH

Dispense the required amount of Prisma APH on a transfer block and protect it from ambient light.

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Fill the cavity with a small surplus, allowing the material to extend beyond the cavity margins.

POSTERIOR TEETH

On the posterior teeth, incremental application and polymerization is recommended in order to minimize shrinkage of the polymerization. For Class I and II cavities polymerize each increment (2mm thick or less) separately.

POLYMERIZATION

Expose each restoration area to the visible light for at least 20 seconds. For Class I, II, III and IV restorations, the composite should receive further polymerization through the walls of the buccal and lingual/palatal enamel.

FINISHING and POLISHING

The finish can be done immediately after polymerization. Thicker excesses can be removed with multi-blade drills and refinement done using Enhance System tips (cups, discs and points). A high shine is obtained with the use of Prisma Gloss and an Enhance polishing cup.

CONSERVATION CARE: Store at room temperature between 15° C and 30° C and away from direct sunlight.

DATE OF MANUFACTURE AND BATCH: see label

SHELF LIFE: 36 months

PROFESSIONAL USE