Whitegold Whitegold  
office protector

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## WHITEGOLD OFFICE & WHITEGOLD PROTECTOR Blue

In-office teeth whitening system

**WHITEGOLD OFFICE** is a powerful whitening hydrogen peroxide 35% based agent indicated for exclusive use in an office, as sole treatment (whitening in an office), or associated with Whitegold Home (10% or 16%) to accelerate the whitening procedure (combined whitening). The self-mixing system, provided by the syringes assures perfect concentration of the peroxide and its efficiency at each application. Once mixed, phases A and B are ready for use, and can be kept in a refrigerator for a period of 20 days.

**WHITEGOLD PROTECTOR Blue** gingival barrier is an opaque blue color resin easy to apply. When light-cured, becomes hard, thus providing protection to the gingival tissue. It is indicated as an indispensable procedure prior to the application of the Whitegold Office whitening gel, and it is a practical and efficient alternative to using absolute isolation. It can be used either in more extensive cases (several teeth) or when just only or more teeth are involved in the procedure. It is recommended to associate its use to a lip retractor and an aspirator.

### PRESENTATION FORMS:

- Kit content 02 Phase A syringe containing 1.6g + 02 Phase B syringe containing 0.4g + 01 Syringe WHITEGOLD PROTECTOR Blue containing 2.4g + 05 Application Tips
- Refill content 01 Phase A syringe containing 1.6g + 01 Phase B syringe containing 0.4g + 02 Application Tips
- Package content 01 WHITEGOLD PROTECTOR Blue syringe containing 2.4g + 02 Application Tips

### COMPOSITION:

- WHITEGOLD OFFICE PHASE A: Hydrogen Peroxide and Avicel RC-591
- WHITEGOLD OFFICE PHASE B: Tickner gel composed by Glycerin, Carbolpol 974 NF, Special mint flavor, Sodium Hydroxide solution 50% and Solution 1% Strain FD&C Green.
- WHITEGOLD PROTECTOR Blue: Methacrylate based resin, BHT, Inert load, Tertiary amine, Stain and Camphorquinone.

### INDICATIONS:

- WHITEGOLD OFFICE:

This is indicated for vital teeth with extrinsic color alteration with chronic character, for which conventional procedures are not efficient (prophylaxis) and for teeth with less severe intrinsic color alterations (physiological darkening by dentine placement), traumas, etc. For more severe intrinsic color alterations, such as fluorosis, severe staining by tetracycline, although effective, and the whitening results are less foreseeable. The use of whitening agents is counter-indicated in cases of imperfect Amelogenesis and dentinogenesis. Whitegold Office can also be indicated as a complementary technique for non-vital teeth whitening, associated to the vitalized teeth whitening.

### - WHITEGOLD PROTECTOR Blue:

This is indicated for EXCLUSIVE use with the vital and non-vital teeth whitening system, through a professional procedure in a dental practice office. For non-vital teeth, the "Walking Bleach" technique is not recommended. The use indication shall be accompanied by isolation for intra-oral protection, and the patient's and the professional's eyes protection (safety goggles).

### CONTRAINDICATIONS:

- WHITEGOLD OFFICE:

Its application is counter indicated on restorative material with the purpose to promote any kind of whitening. It is recommended to carry the whitening treatment first, and replace the restorations in two weeks from its finalization. This is due to the fact that there is a reduction in the adhesive efficiency provided by the peroxide presence and release and the need of color stabilization, which might achieve higher esthetic results.

- WHITEGOLD PROTECTOR Blue:

Avoid applying the product on the gums, in case it is not in perfect health conditions. Ex: in a period of post-surgical intervention healing; in the presence of irritation or gingival inflammation.

### WARNINGS:

- WHITEGOLD OFFICE:

- 1- Whitegold Office shall not be used in patients with history of allergy to Hydrogen Peroxide;
- 2- Its use is not recommended in pregnant or breast-feeding patients;
- 3- The gel contact with the gums or mucosa may cause irritation or injury due to its caustic action;
- 4- The product is designated to topic application, thus it shall not be swallowed;
- 5- The gingival barrier shall not be used in patients reporting allergy to resin. It would recommendable to test on a small area of the gums before its application;
- 6- Whitegold Office is a system for exclusive use in an office, this it shall be applied under professional supervision.

### PRECAUTIONS:

- WHITEGOLD OFFICE:

1. The professional shall be duly equipped whenever working with peroxide (goggles with side protection, a jacket and gloves);
2. The patient shall be protected since the beginning of the procedure till the removal of the product, using goggles with side protection;
3. Be careful when removing the syringe protection tip. During the activation (manipulation) of the product, it shall be kept far from the patient. The same cares shall be taken when connecting the syringes after mixing.

- WHITEGOLD PROTECTOR Blue:

1. Incomplete polymerization of the product or yet excessive time might lead to slight tissue irritation.
2. In case the patient has more severe allergic reaction to the product, it is sufficient to interrupt the use of the product;
3. Non-injectable product. Do not reuse the syringe.

### ADVERSE REACTIONS:

- WHITEGOLD OFFICE:

Some patients may present certain sensitivity, mainly in cases of root exposure, pre-existing restorations, presenting micro-leakage or dentine exposure. Before the beginning of the treatment, the patients have to be warned about the possibility of transitory sensitivity occurrence.

- WHITEGOLD PROTECTOR Blue:

The product contains monomers, which might provoke slight allergic or irritating reactions in patients sensitive to these substances. During the product polymerization, slight feeling of heat or burning on the gums might be experienced as a result from the product polymerization and contraction reaction.

### SPECIAL CARES:

- WHITEGOLD OFFICE:

In case of contact of Whitegold Office with the oral mucosa or on the skin, immediately wash with plenty of water. Apply an ointment of sodium bicarbonate and water on the affected area or in case there is no ointment, prepare a neutralizing solution in a dappen pot or in a flask with a lid, mixing sodium bicarbonate powder with physiological serum till forming a homogeneous liquid phase. It is possible that a sodium bicarbonate precipitate occurs in the prepared solution on the bottom of the vessel after mixing in the form of a paste, which can also be used. If the symptom persists, consult a physician.

### USE INSTRUCTIONS WHITEGOLD OFFICE:

- 1- Removal of calculus with the help of ultrasound Cavitrón associated to the cleaning of the dental surfaces with pumice and water. This enabling either the removal of the plaque or the extrinsic stains, guaranteeing the whitening gel action and efficiency;
- 2- Record the initial color of the teeth before submission to the whitening procedure with the help of the strain scale and record the initial color on a card;
- 3- Isolation of the teeth to be whitened - done with gingival barrier WHITEGOLD PROTECTOR Blue associated to the lip retractor or under absolute isolation. For the application of the gingival barrier WHITEGOLD PROTECTOR Blue, see instructions for use in subsequent topic.
- 4- Application of the proper whitening agent: Before use, keep the gel at room temperature (22°C). Remove the protection lid of the Phase A syringe, and then do the same with the Phase B syringe. Connect both syringes, screw till feeling their final lock. Inject the content of the Phase A syringe inside the Phase B syringe and repeat the same operation in the opposite direction till completing 20 times, i.e., 10 movements in each direction softly. After completing the mixture, return the gel to the Phase A syringe.  
Disconnect the syringes carefully and adapt the application tip to the Phase A syringe. Poor the gel on the dental surfaces with sweeping movements. In case air bubbles are detected during application, remove the syringe from the patient's proximity and eliminate the air bubble on a surface before continuing with the product application. Make sure that all surfaces in question are completely covered with gel;
- 4.1- Whitening of associated non-vital tooth - in case of non-vital teeth associated to vitalized teeth, first it shall be proceeded with crown opening of the devitalized tooth, cleaning the chamber for the removal of the restoring material and/or possible endodontic material remnants, preparation of the cervical "plug", application of the whitening gel Whitegold Office either on the vestibular surface or inside the crown cavity. As in the case of the non-vital teeth, the gel can be kept for up to 45 consecutive minutes with no need to change in the same session. After the gel removal (washing), it shall be proceeded with the cavity temporary sealing.

### TECHNIQUE RECOMMENDED BY THE MANUFACTURER

- 1- The manufacturer suggest to maintain the product in contact with the dental surfaces for a period of 30 to 45 consecutive minutes (upon the professional's criterion), with no need to change consecutively at every 15 minutes. This fact is due to the properties of the product to maintain its pH neutral for a longer period of time after mixing.
- 2- Application of light source is not necessary as an additional method to accelerate the whitening process. The product will act efficiently, including without light activation. In fact, this light source activation, according to many authors, can accelerate the surface whitening by means of breaking extrinsic stains; however, it cannot act on intrinsic stains, because, for this purpose, longer period of time would be needed for the gel diffusion.

- As all whitening procedures, the color change is variable and depends directly on the dental permeability. In addition to producing whitening due to the perfect gel viscosity and the consequent penetration power as it maintains its pH neutral for longer periods, Whitegold Office can be certainly kept for up to 45 minutes with no need of consecutive changes, thus simplifying the procedure.

- In order to obtain better results, 2 to 6 visits average are necessary, depending on the case and on the response to the whitening treatment. In patients with sensitivity history or potential, the time can be reduced to 30 minutes application, increasing the number of visits.

### USE INSTRUCTIONS WHITEGOLD PROTECTOR Blue:

1. First remove the protection lid of the syringe, then connect the application tip;
2. After the prophylaxis of the involved teeth, position a lip and cheek retractor, which is more compatible with the dimensions of the patient's mouth;
3. Wash and dry the teeth well, paying attention for a 4-6mm range of the gums to be completely dry, in order to enable the gingival barrier effective adhesion;  
Remark: When the isolation is in the lower region, the saliva flow control is more critical and shall be done by means of cotton rolls, placed in the sublingual region associated to the aspirator. In this case, the isolation of the region can be complemented by extending the barrier to the cingulum and embrasure region of the front lower teeth, thus, preventing from saliva over-leakage towards the dental surface.
4. With the syringe tip positioned the most parallel possible to the vestibular face of the tooth in question, apply a 1-2mm thick layer along the whole extension of the gums to be protected. Make sure to have covered the whole gingival margin, extending 0.5mm forwards the enamel surface and at least 4-6mm in the apical direction. Perform the final visual inspection, observing whether there are possible failures, which have to be totally filled in;
5. Once the correct application of the product is assured, light cure along the whole extension of the barrier, avoiding the equipment tip to remain in one area only for more than 2 seconds. A barrier with extension involving 8 teeth needs 20 seconds of cure average or less, when 600mW/cm2 light intensity is used. Shorter period will be necessary, if the device intensity is higher than the mentioned.  
Caution: Excess of exposure in a localized area can produce excessive heat and consequently generate tissue injury. We recommend the activation of one quadrant first, and then the activation of the other.
6. Test of the barrier surface hardness/resistance with the help of a dental explorer. When adequately polymerized, the barrier presents shiny aspect. If any gingival area is still exposed, apply an additional layer/till its complete coverage and photopolymerize;
7. For removal after the completion of the whitening session, first make sure that the gel has been completely removed from the dental surface. After that, with the help of a dental explorer raise the end of the barrier and pull with the tweezers, detaching it completely.

### STORAGE AND TRANSPORT:

- WHITEGOLD OFFICE:

Keep the product at 5°C to 20°C temperature, considering that it can be transported at temperature up to 30°C, and far from reach of children.

-WHITEGOLD PROTECTOR Blue:

Keep in a place protected from direct light, heat and far from reach of children. Store the product at 10°C to 25°C.

### BATCH, DATE OF VALIDITY

See the external packaging.

### VALIDITY PERIOD:

- WHITEGOLD OFFICE: 30 MESES
- WHITEGOLD PROTECTOR: 24 MESES

Reprocessing prohibited.

### PROFESSIONAL USE