EVOLU-X - Restaurador Nano-híbrido Fotopolimerizável

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EVOLU-X - Light-Cured Nanohybrid Restorative

The Nanohybrid restorative EVOLU-X® is a radiopaque composite activated by visible light which can be used in anterior and posterior restorations both primary and permanent teeth. Before its application a suitable dental adhesive such as XPBond™, or Prime&Bond® 2.1 must be used. This restorative material comes in syringes with 10 colors of the Vita® scale (T, XL, C2, BW-O, B1, B2, A3.5, A3, A2, A1).

COMPOSITION

EVOLU-X Raw Material	
Silanized barium aluminum boro silicate glass	Dyes
Silanized barium fluoro aluminum boro silicate glass	Photoinitiators
Dimethacrylate Bis-GMA and BIS-EMA	Inhibitors
Nanoparticulated Silica	

1- Indication, purpose or use

The restorative material Evolu-X is indicated for use as a direct restorative of anterior and posterior teeth. It can also be used as restorative material in esthetic cases, such as closing of diastema and lengthening of incisors, transformation of conoid teeth, veneirs and the indirect manufacturing of inlays and onlays.

TECHNIQUES OF USE

MODE OF USE

STEP-BY-STEP INSTRUCTIONS - DIRECT RESTORATION

1. Selection of color

Before selecting the color, the teeth must be clean, hydrated and free of foreign material or stains.

2. Preparation and Insertion

After performing the preparation of the cavity and applying the adhesive system (XPBond or Prime & Bond 2.1). Remove the top of the syringe and dispense the quantity required of the Evolu-X material turning in a clockwise direction. To prevent the material overflowing after use, place the top, point the syringe upwards and turn it in a counterclockwise direction.

3. Polymerization

Insert in the cavity increases of up to 2mm and photoactivate each area of the surface of the restoration with a dental curing light appliance (wavelength of 470 nm). The minimum light power must be about 300mW/cm2 for 20 seconds.

4. Finishing and Polishing

- 4.1. Start the finishing immediately after the process of photoactivating the last layer. Remove the excess material and establish a general outline with finishing burs made of carbide or fine diamond tipped. Furthermore, it is advisable to use the Enhance® finishing systems. See the manufacturer's instructions for use.
- 4.2. In order to achieve a great surface shine on the Evolu-X restorative material, it is necessary to execute the polishing. We recommend the use of the micro-polishing system with diamond PoGo® and/or the polishing pastes such as Prisma Gloss composites. See the manufacturer's instructions for use.

CLEANING AND DISINFECTION

In order to avoid that the syringes with the restorative material Evolu-X are splashed by fluids from contaminated hands or body, or oral tissues, it is advisable to use a protective barrier, to avoid the syringe being contaminated.

2. Precautions, Restrictions, Warnings, Special Care, Storage and Transportation

Warnings

1- The restorative material Evolu-X contains monomers of polymerizable metha crylate. Avoid prolonged or repeated contact with the skin, oral tissue and eyes. Avoid prolonged inhalation. Do not swallow.

Contact with the eyes: the restorative material Evolu-X contains methacrylates which can irritate the eyes. When using this product, use protection for the clothes, eyes and hands of the patient and the professional.
If there is contact with the eyes, wash them with plenty of water and seek medical aid.

Contact with the skin: the restorative material Evolu-X contains polymerizable monomers which can cause irritation of the skin (contact allergic dermatitis) in sensitive people. If there is contact with the skin, clean it immediately with cotton wool soaked in alcohol and then wash the skin with water and soap. If there is irritation or another allergic reaction, stop using the product and seek medical aid.

Contact with oral mucosa: avoid the contact of the non-polymerized material with the oral tissue. If there is accidental contact, wash the affected area with plenty of water. If there is irritation of the mucous, seek medical aid.

Swallowing: do not inhale or swallow. If swallowed by accident, drink plenty of water. If there is nausea or physical indisposition, seek medical aid.

During the handling of chemical products always use individual protective barriers as stipulated by ANVISA (safety goggles, mask, jacket and cap), as well as for the patient (safety goggles, apron and cap).

PRECAUTIONS

- 1- This product is manufactured to be used only as per the Instructions for Use. The use in a way differing from that described in the instructions for use is solely and exclusively the responsibility of the professional who does so.
- 2- The contact with the saliva and blood during the application of the composite can harm the restoration. It is advisable to use absolute isolation, cotton wool or rubber protection or another suitable technique of isolation.
- 3- Wear a jacket, mask, gloves and safety goggles for appropriate protection. It is advisable that the patients wear safety goggles.
- 4- The restorative material Evolu-X is a material which polymerizes with light. Once out of the packaging, it should be used quickly or else protected from ambient light.
- 5- The material should come out of the syringe easily, exerting a slight and Constant pressure. DO NOT APPLY TOO MUCH FORCE.
Excessive pressure can make the material come out in an unforeseen manner.
- 6- The syringes should be hermetically closed after use.

RESTRICTIONS OR CONTRAINDICATIONS

Reactions to methacrylate resins.
All the clarification required and care with use are in the instructions for use for the product.

STORAGE AND TRANSPORTATION

Store at a temperature between 10°C and 24°C. Protect from humidity. Do not freeze. Do not use after the validity date has expired.

ADVERSE EFFECTS

The non-polymerized material can irritate the eyes and skin.

- Contact with the eyes: Irritation and possible damage to the cornea.
- Contact with the skin: Irritation and possible allergic reaction. Red stains may occur on the skin.
- Mucous: Inflammation and scaling.

INTERACTIONS WITH OTHER MATERIALS

The material should not be used with dental materials which contain eugenol, as they could interfere with the polymerization process and cause a reaction of the components of the material.

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3. Manner of presenting the product

Syringe containing 4g, in the colors T, XL, C2, B1, B2, BW-O, A3.5, A3, A2, A1.

Expiration Date

3 years. Do not use after expiration date.